



ENDURO POWERLITE ^{IK2PCS}

FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

CARBON FIBRE

12.3kg

**PRODUCT
WEIGHT**
(without battery)

**ON BOARD
CONTROLLER**



**ANTI-TIP
WHEELS**

FEATURES

- Carbon fibre frame
- Folds up for transport and storage
- Comfortable padded seat and backrest
- Fold down backrest with storage pocket
- Lift up armrests
- Onboard controller with LCD display and fingertip controls
- Battery power indicator
- Variable speed setting
- Rear wheel drive
- Front wheel suspension
- Electromagnetic brake
- Airplane safe battery
- Rechargeable Lithium-ion battery
- Anti-tip wheels
- Calf strap
- Lap belt
- Horn



**AIRLINE
TRAVEL
READY**

ASK FOR THE BRAND

ikonmobility.com.au



ENDURO POWERLITE IK2PCS

FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

CARBON FIBRE



**MOTORISED OR
FREEWHEEL**



**FOLDS UP
FOR STORAGE**



**CONVENIENT
STORAGE BAG**

**FLIP UP
ARMRESTS**



**REMOTE
CONTROL
INCLUDED¹**

SPECIFICATIONS

- Overall (L to unfolded footrest): **1045mm(L) x 555mm(W) x 885mm(H)**
- Folded (L to folded footrest): **820mm(L) x 275mm(w) x 800mm(H)**
- Seat (H from floor): **420mm(L) x 460mm(W) x 490mm(H)**
- Backrest (H from seat): **440mm(W) x 420mm (H)**
- Push Handle Height: **880mm**
- Turning Circle: **1080mm**
- Speed: **6km/hr**
- Max Range: **Up to 18km²**
- Safe Climb Angle: **6°**
- Max Climb Angle: **6°**
- Motor: **2 x 250W high torque brushless motors**
- Battery: **24V / 10.4Ah Lithium-Ion**
- Charger: **2 Amp Lithium-Ion**
- Tyres (flat free): **Front: 165mm (6.5") / Rear: 280mm (11")**
- Product Weight: **13.3kg with battery, 12.3kg without**
- Max User Weight: **140kg (on flat ground)**

¹ Limited time offer.

² Dependant on occupant and conditions.

ASK FOR THE BRAND

ikonmobility.com.au