



# Novis Avalon Rotating Bed



220 kg



3  
ELECTRICAL  
WARRANTY  
FRAME  
10



2  
COVER  
WARRANTY  
BODY  
6

The Avalon Rotating Bed is an electrically powered bed and rotating profiling chair suited to any care environment.

The bed offers greater independence to weight bearing patients who require assistance transferring between standing and sitting.

1

Allows user to get out of bed with minimal to no assistance

---



2

Dual function - fully functioning electrical bed that can also be used as a comfortable day chair

---

The Avalon Rotating Bed offers greater independence and mobility to users who are weight bearing but require assistance to transfer between sitting and standing.



3

Removable arm rests to aid transfer in and out of bed

---

5

Backrest auto regression function minimises risk of pressure injury

---



6

Flexible design allows the bed to be configured to rotate to either the left or the right hand side

---

4

Single button control allows easy operation to reposition patient from lying down to seated position

---



## Specifications

CODE	BE1007
SAFE WORKING LOAD	220 kg
WARRANTY	Bed - 10 year bed frame/ 3 year electrical/ 1 year castor Mattress - 6 year body/ 2 year cover
FRAME MATERIAL	Powder coated steel
OVERALL DIMENSIONS (MM)	2324 (L) x 1050 (W)
HEIGHT RANGE (MM)*	500 to 900 mm
MATTRESS	Custom two piece castellated foam mattress Multi stretch, vapour permeable cover 360° opening zipper Adjustable straps for securing to bed frame
CASTORS	100mm lockable twin castors
MATTRESS SIZE RECOMMENDED (MM)	<b>Main Body Section</b> 1750 (L) x 880 (W) x 150 (H) maximum <b>Foot Section</b> 300 (L) x 880 (W) x 150 (H) maximum
BED WEIGHT INCLUDING SIDE RAILS AND ARM REST	149 kg
BACKREST ANGLE	0 to 70°
KNEE BREAK ANGLE	0 to 35°
TRENDELENBURG TILT	0 to 11°
ANTI-TRENDELENBURG TILT	0 to 11°
COMPLIANCE	IEC 60601-1:2005 IEC 60601-2-52:2012 IEC 61000-3:2008 IEC 61000-4:2009
ARTG	Bed - 278782 Mattress - 176837

\*measured from top of mattress platform

## Accessories

BE7003	Self Help Pole
BE9006-WH	Fitted Sheet Set, White

