

Shoprider

778SL

ASK FOR THE BRAND **SHOPRIDER**



SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

778SL

SHOPRIDER[®]
www.shoprider.com.au *mobility*

The SHOPRIDER™ 778SL is a versatile performer from the impressive SHOPRIDER™ range of scooters. The 778SL is loaded with great features to provide the user with outstanding functionality and performance including a powerful motor, large battery capacity, suspension, rotating / sliding seat, modesty panel with dual driving lights & indicators. The reliability of this scooter is further enhanced by the rear wheel drive and electromagnetic brakes.

SPECIFICATIONS

778SL

Length	1230mm
Width	590mm
Height	1250mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	136Kgs
Charger	4Amp
Motor Size	1.2Hp
Maximum Speed	12Kph
Battery Capacity (Minimum Amp)	2 x 50Amp
Maximum Range #	up to 30km
Wheel Size	260/330mm
Turning Radius	1250mm
Weight (w/o batteries)	76kg

FEATURES

High Back Seat
Adjustable Steering Tiller
Adjustable Width Arm Rests
Modesty Panel
Fold Arm Rests
Sliding Rotating Seats
Front Shopping Basket
Variable Seat height
Battery Level Warning Light
Front & Rear Suspension
Front Lights
Anti Tips Wheels
Comprehensive Warranty

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia and is manufactured to comply with all relevant international ISO standards.


IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING: Do not exceed 10kph on public roads or paths

#Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

778SL



SHOPRIDER[®]
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au