

VICTORY[®] 10 LX WITH CTS SUSPENSION

4-WHEEL SCOOTER

THE ULTIMATE MID-SIZE LUXURY SCOOTER

The **Victory[®] 10 LX with CTS Suspension** Mobility Scooter features CTS (Comfort-Trac Suspension) Suspension and brings a new level of performance and features to mid-size scooters. Featuring CTS independent front and rear suspension, sleek style, feather-touch disassembly and standard lighting, the 4-wheel Victory 10 LX with CTS Suspension is the ultimate mid-size luxury scooter.



Front Suspension



Rear Suspension



PRIDE[®] MOBILITY SCOOTERS
LIVE YOUR BEST[®]

www.pridemobility.com.au

SPECIFICATIONS

VICTORY® 10 LX WITH CTS SUSPENSION

4-wheel

Weight capacity:	181.4 kg (400 lbs.)
Maximum speed ¹ :	Up to 9 km/h (5.6 mph)
Ground clearance ² :	5.6 cm (2.2") at motor
Turning radius ² :	152.1 cm (59.9")
Overall length ² :	119.4 cm (47")
Overall width ² :	56.5 cm (22.2")
Seat-to-floor height:	45.1 cm (17.7")
Front tyres ² :	9 cm x 25.4 cm (3.5" x 10") pneumatic; 25.4 cm (10") pneumatic
Rear tyres ² :	9 cm x 25.4 cm (3.5" x 10") pneumatic; 25.4 cm (10") pneumatic
Suspension:	Front and rear
Range per charge ^{1,4} :	Up to 26 km (16.16 miles)
Maximum allowable incline:	6°
Total weight without batteries ⁶ :	70.9 kg (156.2 lbs.)
Heaviest piece when disassembled:	35.2 kg (77.6 lbs.) (front section)
Standard seating:	Type: Anskon Junior seat Weight: 19 kg (42 lbs.) Material: Black Vinyl Dimensions: 43.16 cm (17") width, 43.16 cm (17") depth
Drivetrain:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes:	Regenerative and electromechanical
Battery requirements ⁵ :	(2) 12-volt, deep cycle Size: 40 AH (standard) Weight ³ : (40 AH) 9.7 kg (21.4 lbs.) each
Battery charger:	5 amp, off-board charger
Color:	  Candy Apple Red Viper Blue

FEATURES

- CTS Suspension (Comfort-Trac front and rear independent suspension)
- Sleek, modern design
- Exclusive black, non-scuffing tyres
- Feather-touch disassembly permits simple frame separation with only one hand
- Modular design for easy serviceability
- Delta tiller with wraparound handles
- Easily accessible charger port in tiller
- Full LED lighting package that includes curb and rear lights and directional signals
- Courtesy light under tiller
- Storage compartment below tiller
- Enhanced flush-mounted horn buttons
- Auto-connecting front-to-rear cable

OPTIONS

- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Saddlebag
- Rearview mirror

- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).
- ² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ³ AGM or gel cell type required.
- ⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁵ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.

AU_VICTORY10LXCTS-SS-08-17-15

Pride
Mobility Products
Australia Pty. Ltd.
— LIVE YOUR BEST —

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®



(Tel) 03 8770 9600 • (Fax) 03 9703 2720
20-24 Apollo Drive Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 661

www.pridemobility.com.au

