


Etac Contour pen and Uni turner

We all need to write small notes or sign a document from time to time. There are also knobs and handles that need to be turned in order to open doors, cupboards or switch lamps on/off.

Etac Contour pen


The Contour pen is specially designed to enable handwriting with poor gripping ability. Thanks to the narrowed grip area the pen rests comfortably in the hollow of the thumb and index finger. The thumb stabilizes the grip and writing is possible with very little effort. The light fluid ink further reduces the strain on the fingers.



The Contour pen lays comfortably in the hollow of the thumb and index finger.

Etac Uni turner

Uni Turner makes small, day-to-day tasks easier, such as opening and turning small knobs and handles. It has an ergonomic and comfortable handle which facilitates pressing the knob before turning. The red hexagonal pegs ensure a good grip around the knob.



The Uni Turner offers a firm and comfortable grip for the wrist and relieves the finger joints when in use.



Description

Contour pen
Uni universal turner

Item no.

80701001
80503003

Weight

Contour Pen: 26 g (1 oz)
Uni turner: 84 g (3 oz)

Material

Contour pen: Handle ABS
Uni turner: ABS
Pegs: Polyamide

Design

Contour Pen: Designkonsulterna
Uni turner: Ergonomidesign