

# MAVERICK 12



**ScooterWorld**  
AUSTRALIA

**merits**  
*with you for life*

[www.scooterworldaustralia.com.au](http://www.scooterworldaustralia.com.au)

# MAVERICK 12



## FEATURES & BENEFITS

- Mid-wheel drive design delivers outstanding maneuverability in tight, compact spaces and a smooth and stable ride both indoors and outdoors, even over tough terrain
- Patented obstacle climbing suspension that works
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Low centre of gravity design delivers a smooth and stable ride
- Swing away joy stick
- Height and width adjustable armrests
- Height adjustable four post seat with reclining back rest

|                           |  |
|---------------------------|--|
| Range                     | 25 - 35km  |
| Motor                     | DC 24 Volt / 575W                                |
| Controller                | Dynamic Shark 50A                                |
| Battery                   | 12V/40-55 a/h x2 pieces                          |
| Charger                   | 24V/5 amp off-board                              |
| Gradient                  | 8°   |
| Caster Wheel              | F:6" R:6" foam filled tyre                       |
| Drive Wheel               | 12" pneumatic tyre                               |
| Weight capacity           | 205kg  |
| Brake                     | Intelligent, regenerative, electromagnetic brake |
| Length (Footplate down)   | 90 cm  |
| Width                     | 69 cm  |
| Seat Width                | 51cm   |
| Seat Height (from deck)   | 43-48cm  |
| Seat Height (from ground) | 56-58cm  |
| Turning Radius            | 53cm   |