

Back-Eze

Back Support Cushions



The Back Eze is a lightweight, polyurethane moulded foam back support cushion, contoured both vertically and horizontally to support the spine and lumbar area.

The vertical curve helps the spine to maintain the correct “S” shape (lumbar lordosis). The horizontal curve gives firm, gentle support against “rolling” of the body left or right.

Designed for use in the car, office or wheelchair the Back Eze is available in two models – **Standard** and **Mark II** – with differences in the amount of contouring. The Standard model suits most people and is recommended for use in a car. With deeper sides that hug the body and provide added lateral support, the Mark II is designed for a larger user and also suits people with a smaller curve in their lower back.

Includes a navy blue stretch fabric cover that is easy to sponge clean. May be used whether sitting or lying. For best results it should be used against a fairly flat, firm surface.

ITEM	Back Eze Standard	Back Eze Mark II	Replacement Cover for Mark II
CODE	Standard-Back-Eze	Mark-II-Back-Eze	Back-Eze- Replacement-Covers
WIDTH	350 mm	350 mm	
HEIGHT	320 mm	330 mm	
DEPTH	70 mm	110 mm	
COLOUR	Blue	Blue	
WEIGHT	450 g	500 g	
WARRANTY	1 year		

