

Introducing ProCair
by Novis – alternating
solutions for the prevention and
treatment of pressure injuries and other
preventable complications of immobility.

	Standard	King Single	
CODE	PTAM8210	PTAM8211	
CAPACITY (SAFE WORKING LOAD)	40 - 250 kg		
SYSTEM ALTERNATION	1 in 3		
WARRANTY	3 year		
CARRY BAG	Deluxe carry bag with padded pump and mattress compartments		
PRESSURE INJURY RISK	High to very high risk		
COMPLIANCE	IEC 60601-1, IEC60601-1-2, IEC60601-1-11		
ARTG	176273		
CARTON DIMENSIONS (L x W x H) in mm	1125 x 430 x 370 mm		
DIMENSIONS (L x W x H) in mm	2000 x 885 x 200 mm	2000 x 1050 x 200 mm	
MATTRESS WEIGHT	11.5 kg	12.5 kg	
CPR TYPE	Pull-tag type, located at head end		
CELLS	NUMBER	19	
	STATIC HEAD CELLS / ZONE	3	
	NARROW FOOT CELLS / ZONE	5	
	CELL-IN-CELL	Yes	
	CELL MATERIAL	Pure TPU	
	MATTRESS AIR SUPPLY LINE	Flexible silicone piping	
TOP COVER	MATERIAL / TYPE	Blue, polyurethane-laminated nylon – water resistant, vapour permeable multi-stretch	
	FIRE RETARDANCY	BS7175 Crib 5	
	MAXIMUM WASH TEMPERATURE	95° C	
	ATTACHMENT	360° zipper	
	SEAM CONSTRUCTION	RF welded on all external seams	
BASE	MATERIAL / TYPE	Black PU-laminated nylon, fully sealed	
	HINGE	Yes	
	FIXING STRAPS	4 x webbing at foot end and hinge area	

Control
Unit

DIMENSIONS (L x W x H) in mm	120 x 315 x 270 mm	
WEIGHT	4.5 kg	
POWER	240V 50 Hz	
COMPRESSORS	Dual	
OPERATION	Full digital control	
CYCLE TIME	12 minutes	
NOISE LEVEL, WITHIN 2 METRES	38dB(A) ± 2	
MIN / MAX PRESSURE	18-60mm HG +/- 6mm HG	
ELECTRICAL PROTECTION TYPE	Class II Type BF	
OPERATING ENVIRONMENT	Operation 15°C to 35°C Storage 5°C to 60°C	
PRESSURE SETTING	AUTOMATIC	NA
	MANUAL	8 levels
FOWLER FUNCTION	No	
BACKUP BATTERY	No	
ALARMS	Low Pressure High Pressure Alternation Failure Startup Failure Power Failure	



THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.