

Shoprider

RAINRIDER

ASK FOR THE BRAND SHOPRIDER



SHOPRIDER[®]
www.shoprider.com.au
mobility

RELIABLE MOBILITY GUARANTEED

Shoprider

RAINRIDER

The SHOPRIDER™ RAINRIDER is the only fully covered mobility scooter currently available in the world. The RAINRIDER is perfect for person who needs to be protected from the bad weather & in summer the doors easily remove for the operators comfort. The RAINRIDER comes with 330mm Wheels, Double Pillow Comfort Seat, 2hp High torque motor, Lockable rear parcel carrier and all the features you expect from the SHOPRIDER range, rotating lockable seat, suspension and adjustable steering tiller.

The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

	<u>RAINRIDER</u>
Length	1600mm
Width	790mm
Height	1620mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
*Maximum Load	Single user 160Kgs
Charger	8Amp
Motor Size	2.0Hp
Maximum Speed	10Kph
Battery capacity	up to 2 x 75 Amp
Maximum range#	up to 40km
Wheel Size	330mm
Turning Radius	1800 mm
Weight (with battery)	195kg

WARNING: Do not exceed 10kph on public roads or Paths

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.


Driving range varies with the factors shown below

- (1)The weight of the occupant
- (2)Level ground & Ambient temperature
- (3)Battery condition / Type of charger
- (4)Driving Style

FEATURES

Double Pillow Comfort Seat
Delta Steering Tiller
Adjustable Steering Tiller
Lockable rear parcel carrier
Sliding Rotating Seat
Variable Seat Height
Battery Power indicator
Double Head Lights
Rear Lights & Indicators
Hazard Lights
Windscreen wiper
High & Low Speed
Removable doors
Comprehensive Warranty

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)



12" (20) – 70kg
9" (15) – 90kg
6" (10) – 110kg
0" – 160kg



For your nearest Distributor call 1800 672 633
or visit our Web Site at www.shoprider.com.au

Your Authorised dealer