

Scooter S941A

4-Wheel Full Suspension Electric Scooter

- Brand new sporty streamline design
- Full suspension system with shock absorbers delivers a smooth and comfortable ride
- Modern LED lighting system delivers a brilliant light with less power consumption
- Digital LCD dashboard display: time display, battery gauge, temperature, maintenance messages, speed display and odometer
- Extra storage space beneath tiller
- Automatic safety turning half speed control by sensor




LED lighting system



LED lighting system



Digital LCD dashboard display





More space for storage

Technical Specs.

Weight Capacity	159kg/ 350lb
Max. Speed	12kph/ 7.5mph; 15kph/ 9.5mph
Range	50km/ 32mi.
Turn Radius	167cm/ 65.7"
Ground Clearance	9cm/ 3.5"
Motor	DC 24V, 400W(Rating)/ 1750W(Max)
Controller	120A
Battery	12V/ 75AH x 2pcs
Charger	5A off-board
Gradient	up to 10°
Front Wheel	12"
Rear Wheel	12"
Brake	intelligent, regenerative, electromagnetic brakes

Dimensions



seat dimensions [cm/in]							
D	E	F	G	H	I		
46/18	40/16	46/18	63.5/25	43-48/17-19	97/38		
overall dimensions [cm/in]			weight [kg/lb]				
A	B	C	base	battery/ea	seat	ttl w/batt.	capacity
130/51	69/27	125/49	69/152	25/55	18/40	118/260	159/350