

## OVERBED TABLE

- Quality laminated timber top.
- Protection side strip.
- Multi angle adjustable tabletop.
- Compact folding for storage.
- Unique tilt action design brings the table top closer to the patient.

IPA468

Laminated timber top

Protection side strips

Angle adjustment knob

Sturdy steel frame  
- powder coated finish

Height adjustment knob

Lockable castors

Tilt adjustment knob

TSA468-B - additional table top



Optional  
Additional  
table top

Table in tilt position



Table top  
- angle adjustable

Table tilts closer towards  
the seated patient

Height adjustable



## PRODUCT PACKAGING



IPA468

overall height	overall length	overall width	castor size	unit weight
min 66cm max 98cm	60cm	48cm	5cm	8.3kg

**W**  
1 YEAR WARRANTY  
\* conditions apply