

Adjustable Seat Height Walker



RG4210ADJ



RG4209ADJ

Features

- Easy multiple seat height adjustment
- Available in 6" or 8" wheel size
- Aluminium powder coated frame
- Available in Red, Blue or Green marble finish colours
- Multi height adjustment handles
- Large soft under seat storage bag
- Comfortable folding backrest
- Individually boxed for easy transport

Dimensions

RG4209ADJ 6" Wheel	
Seat height adjustment range	470 - 560mm
Handle height adjustment range	700 - 900mm
Max user capacity	120kg
Walker net weight	6.5kg
RG4210ADJ 8" Wheel	
Seat height adjustment range	520 - 580mm
Handle height adjustment range	750 - 950mm
Max user capacity	120kg
Walker net weight	7.5kg

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@redgumbrand.com.au