

New Product Release

Space Saver Shower Stool



SPACE SAVER SHOWER STOOL  Australian Design

BATHROOM SAFETY

-  **TOOL FREE**
Easy Assembly
-  **WEIGHT CAPACITY**
130 kg
-  **TESTED AND CERTIFIED**
TO AUSTRALIAN STANDARDS
-  **TGA REGISTERED**

REDGUM Daily Living and Lifestyle Products

Patent Pending Model : RG 554H QTY 1 www.redgumbrand.com.au

Features Include

- Australian Standards Approved
- Polypropylene Reinforced Seat Top
- Patented Leg Locking System
- Adjustable Aluminium Legs with Sure Lock Adjustment & Flexible PVC Leg Tips
- 130Kg Maximum User Weight
- Strong Lightweight Construction
- Lifetime Warranty Against Corrosion
- Easy Assembly (No Tools Required)
- Angled Seat for Maximum Comfort
- Easy Clean Design
- Packaged for Display

Product Dimensions

Seat Width	460mm
Seat Depth	300mm
Seat Handle Height	190mm
Seat Height	490mm to 620mm
Maximum User Weight	130kg
Carton Dimensions	230 x 520 x 310mm
Part Number	RG554H

New products are currently in the processes of being developed using the unique patented leg fitting system

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@shoprider.com.au