

# Afiscooter S3

The Afiscooter-S is winner of the prestigious Red Dot Design Award.

It offers an impressive driving range, an ergonomic design for a noticeably more comfortable ride, a programmable electronic controller, LCD control screen and LED headlights.

Furthermore, the Afiscooter S has a fully adjustable tiller and an integrated canopy for protection from sun and light rain; as well as a fully adjustable orthopedic seat, full suspension and a shock absorbing steel frame.

The three wheel model has superb maneuverability and unrivaled ease of steering.

Available Colours



## **Afikim Electric Vehicles**

has been developing and manufacturing advanced, electrically-powered mobility scooters since 1978! Afiscooters' mobility scooters are designed with reliability and driver safety in mind first and foremost alongside beautiful, and thoughtful design. Afikim Electric Vehicles' worldwide distribution network offers year-round local technical support and meets European ISO 9001-2000, EN 12184 certification, as well as US safety standards with FDA approval.



**AFIKIM**  
ELECTRIC VEHICLES

## Features

- ▶ Fully adjustable swivel orthopedic seat with a movable headrest and armrests
- ▶ Adjustable reach configuration using a single button for forward/reverse
- ▶ Bright, easy to read LCD display control panel indicates speed, battery status, lighting and signaling
- ▶ Powerful, long-lasting, energy efficient LED headlights provide more light and a greater driving range



The sun/rain canopy is winner of the 2011 Red Dot international design award



Improved, secure lockable rear storage  
Front storage for items such as glasses and mobile phone



## Specifications

		AfiScooter-S3 - Single		AfiScooter-S3 - Wide seat	
<b>Overall length</b>		61"	155 cm	61"	155 cm
<b>Overall width</b>		27"	68 cm	--	--
	With wide wheels	35"	88 cm	35"	88 cm
<b>Overall height</b>	Without canopy	49"	125 cm	49"	125 cm
	With canopy	65"	164 cm	65"	164 cm
<b>Ground clearance</b>		5"	12.5 cm	5"	12.5 cm
<b>Weight</b>	With batteries	337 lb	153 kg	359 lb	163 kg
	Without batteries	234 lb	106 kg	256 lb	116 kg
	Add for canopy	+ 26 lb	+ 12 kg	+ 26 lb	+ 12 kg
	Add for wide wheels	+ 22 lb	+ 10 kg	+ 22 lb	+ 10 kg
<b>Travel range</b>		28 m	45 km	28 m	45 km
<b>Weight capacity Inc. driver</b>		450 lb	200 kg	450 lb	200 kg
<b>Maximum speed</b>		9.3 mph	15 Km/h	9.3 mph	15 Km/h
<b>Maximun climbing angle</b>		20%		20%	
<b>Turning radius</b>		56"	143 cm	56"	143 cm
<b>Seating</b>		18"	46 cm	33"	85 cm
<b>Tires type</b>		Pneumatic		Pneumatic	
<b>Tires size</b>	Front	16X7.5-8.0"		3X10"	
	Rear	3.5X10"		16X7.5-8.0"	
<b>Battery type</b>		Two 12V 73A/h		Two 12V 73A/h	
<b>Power unit</b>		Rear wheel drive, 24V DC motor		Rear wheel drive, 24V DC motor	
<b>Color</b>		Metallic Silver			

Tested by manufacturer, based on actual parameters.



**AFIKIM**  
ELECTRIC VEHICLES

**1-800-809-3010** (US only)  
[www.afiscooters.com](http://www.afiscooters.com)

Scan  
the code  
to watch  
the video >>

