

SE

The classic SE is a design milestone for Afikim Electric Mobility and has become an icon in the world of personal electric mobility.

It provides a stable and safe ride in any environment and on any terrain.

The strong steel chassis and rear power unit combine strength with a surprisingly long driving range, as well as an impressive load capacity.

An orthopedic seat offers maximum comfort, and the simple and reliable electronic command system make the SE Scooter extremely easy to operate.

Afikim Electric Vehicles

has been developing and manufacturing advanced, electrically-powered mobility scooters since 1978! Afiscoters' mobility scooters are designed with reliability and driver safety in mind first and foremost alongside beautiful, and thoughtful design. Afikim Electric Vehicles' worldwide distribution network offers year-round local technical support and meets European ISO 9001-2000, EN 12184 certification, as well as US safety standards with FDA approval.



AFIKIM
ELECTRIC VEHICLES

Features

- ▶ A balance of simplicity, reliability and classic design make the SE the world's leading personal electric vehicle for over 30 years
- ▶ Swiveling orthopedic chair for comfort and ease of access
- ▶ Easy to maneuver, quiet and easy operation
- ▶ Choice of selection: single or wide seat, with or without canopy



Extremely sturdy large rear basket



Choice of selection: single or wide seat



Specifications

		SE - Single		SE - Wide seat	
Overall length		62"	158 cm	62"	158 cm
Overall width		29"	73 cm	36"	91 cm
Overall height	Without canopy	49"	124 cm	40"	101 cm
	With canopy	--	--	70"	177 cm
Ground clearance		5"	12.5 cm	5"	12.5 cm
Weight	With batteries	322 lb	146 kg	359 lb	163 kg
	Without batteries	218 lb	99 kg	256 lb	116 kg
Travel range		28 m	45 km	28 m	45 km
Weight capacity Inc. driver		500 lb*	200 kg	500 lb*	200 kg
Maximum speed		9.3 mph	15 Km/h	7.5 mph	12 Km/h
Maximun climbing angle		21%		21%	
Turning radius		45"	115 cm	45"	115 cm
Seating		18"	46 cm	33"	85 cm
Tires type		Pneumatic		Pneumatic	
Tires size	Front	2.75 X10"		2.75 X10"	
	Rear	2.75 X14"		2.75 X14"	
Battery type		Two 12 V 73A/h		Two 12 V 73A/h	
Power unit		Rear wheel drive, 24V DC motor		Rear wheel drive, 24V DC motor	
Rear platform size		--		--	
Color		black & chrome		black & chrome	

Tested by manufacturer, based on actual parameters.

* US Model



AFIKIM
ELECTRIC VEHICLES

1-800-809-3010 (US only)
www.afiscooters.com

Scan the code to watch the video >>

