

Aluminium Shower Stool Clip-On Seat



Suited to smaller shower areas, this mobile shower stool features a 410 mm wide clip-on plastic seat, that is easy to remove for cleaning or storage.

The lightweight frame is corrosion resistant and includes integrated hand grips for increased security and ease of transfers.

The height adjustable legs include rubber tips for added grip; and each leg adjusts independently to match uneven surface requirements.



CODE	BL1140
CAPACITY	110 kg
WEIGHT	2.6 kg
OVERALL WIDTH	530 > 550 mm
SEAT WIDTH	410 mm
OVERALL DEPTH	450 > 490 mm
SEAT DEPTH	365 mm
OVERALL HEIGHT	600 > 730 mm
SEAT HEIGHT	480 > 590 mm
CONFORMS TO	ARTG176831 AS/NZS 3973:2009, AS/NZS 4688.2:2000, AS/NZS 4688.3:2000
WARRANTY	1 year