

Etac Hi-Loo with Brackets

Stable and comfortable toilet seat raiser with brackets.

Offers the best functionality and comfort

Hi-Loo with brackets is a raised toilet seat which is extremely stable and yet easy to remove. The seat is ergonomically shaped and the recesses at the back and front provide for convenient intimate hygiene. Hi-Loo is available both with and without lid and in two different heights.



Etac Hi-Loo with 10 cm (4") and 6 cm (2½") raised seat.

Features



Secure fitting

The brackets are covered with friction rubber and adjustable to fit the shape of the toilet bowl. Hi-Loo fits all common types of toilets.



Easy to lift off

Hi-Loo is easy to fit onto the toilet bowl and can be removed without altering the setting of the brackets.



With lid

Hi-Loo is also available with lid.

Description

Hi-Loo with Brackets, 6 cm (2½")	80301065
Hi-Loo with Brackets, 10 cm (4")	80301105
Hi-Loo with Brackets and lid, 6 cm (2½")	80301067
Hi-Loo with Brackets and lid, 10 cm (4")	80301101

Dimensions

Outer dimensions: depth/width 39 cm/36 cm (15½"/14")
Opening: depth/width 27 cm/20 cm (10½"/7½")

Materials

Seat unit: polypropylene
Brackets: ABS/PC

Weight

6 cm: 1.1 kg (2.4 lbs)
10 cm: 1.5 kg (3.3 lbs)
With lid, 6 cm: 1.9 kg (4.2 lbs)
With lid, 10 cm: 2.3 kg (5.1 lbs)



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.



Etac Hi-Loo Fixed

Stable and comfortable fixed toilet seat raiser.

Etac Hi-Loo with fixed mounting is a discreet toilet seat raiser, that fits most toilets and attaches firmly and securely to the toilet. The ergonomically designed seat has recesses at the back and the front, providing convenient intimate hygiene.

The arm supports are foldable but can be permanently fixed in the folded down position. Hi-Loo can be used with only one arm support if space is limited.

Hi-Loo fixed is available in 6 or 10 cm and with or without arm supports.




Etac Hi-Loo fixed fits most toilets and is easy to clean.

Features



Easy to detach

The snap locks makes it easy to detach the seat unit for cleaning.



Hi-Loo Angled

Recommended when knee or hip mobility is limited.

Accessories

See "Accessories for Hi-Loo and Supporter" on page 50



Description

Hi-Loo fixed, 6 cm (2½")	80301106
Hi-Loo fixed, 10 cm (4")	80301107
Hi-Loo Fixed with arm supports, 6 cm (2½")	80301316-2
Hi-Loo Fixed with arm supports, 10 cm (4")	80301317-2
Hi-Loo Fixed with arm supports, angled	80301318-2

Dimensions

Arm support height from top of seat: 23 cm (9")
 Width between arm supports: 49 cm (19¼")
 Overall width incl. arm supports: 60 cm (23½")
 Hi-Loo Angled rear height: 10 cm (4")
 Hi-Loo Angled front height 6 cm (2½").

Material

Seat unit: polypropylene
 Arm supports: Polypropylene, powder coated steel tubing
 Mountings: aluminium

Weight

Hi-Loo Fixed with arm supports, 6 cm: 4.9 kg (10.8 lbs)
 Hi-Loo Fixed with arm supports, 10 cm: 5.8 kg (12.8 lbs)
 Hi-Loo Fixed, 6 cm: 2.4 kg (5.3 lbs)
 Hi-Loo Fixed, 10 cm: 3 kg (6.6 lbs)
 Hi-Loo fixed with arm supports, angled: 5.5 kg (12.1 lbs)



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.