

Etac OptimaL toilet arm support

Functional and stylish arm supports.

Function and simplicity in a smart combination

Etac OptimaL is a neat wall mounted arm support. The leg ensures stability and the small wall attachment allows fitting close to the toilet. OptimaL stays safely folded and is easily flipped down with the handle.



Accessory



Toilet paper holder
83030001



Description

OptimaL (one piece) 80303006

Dimensions

Length arm supports: 73.5 cm (29")
Floor to arm supports: 80 cm (31½")
Floor to arm supports, folded up: 149 cm (58¾")
Wall bracket: 10.4 x 11.0 cm (4 x 4¼")

Material

Wall bracket, arm supports tubing and supporting leg:
Powder-painted steel tubing
Details: Polypropylene

Weight

3.1 kg (6.8 lbs)



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.