

# Etac Swift Freestanding toilet seat raiser

Lightweight and easy to handle.

## Excellent stability

Etac Swift Freestanding toilet seat raiser provides stable support. Its large footprint and soft, specially designed ferrules give an excellent grip to the floor.

## Functionality

The seat is height adjustable and can be set with a slight anterior tilt. This facilitates standing up and relieves pressure on the hip joints. The legs can easily be taken off for car transportation.

## Flexibility

Besides using the Swift Freestanding over the toilet, it can be used as a shower stool. Its light weight makes it easy to move. With complimentary accessories, the Swift Freestanding may also be used as a commode.



*Etac Swift Freestanding has an integrated splash guard and smooth surfaces which enable easy cleaning.*

## Accessories



### Back support

84005069

Material: Polypropylene



### Soap dish/Lid holder

80209266

Material: Polypropylene



### Soft Seat

84005064

Material: Polyethylene



### Soft Back

84005065

Material: Polyethylene



### Pan

without lid 84005066

with lid 84005068

Material: Polypropylene

## Description

Swift Freestanding toilet seat raiser 81702020

## Dimensions

Total width: 56 cm (22")

Seat width: 54 cm (21¼")

Seat height: 42–57 cm (16½"–22½")

Width between the arm supports: 45 cm (17¾")

## Material

Seat and arm supports: polypropylene

Legs: Aluminum and polyamide

Ferrules: TPE

## Weight

3.9 kg (8.6 lbs)



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.