

Deluxe Transfer Bench

RG9095



The REDGUM Brand Deluxe Transfer Bench has been designed to withstand the rigours of the bathroom environment with its anodised frame, stainless steel fittings and plastic moulded seat. The transfer bench has height adjustable legs and the backrest can be fitted on either side depending on the user's requirements. Two of the legs are fitted with suction cup stoppers for added stability.

Dimensions

- Seat width 800mm
- Leg height adjustment range 450mm to 570mm
- Max user capacity 130kg

Features

- Polyurethane self draining seat
- Removable Backrest that can be fitted either side
- Rust free construction
- Built in soap holder

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@shoprider.com.au