

Baricare Range Aluminium Shower Stool Knock Down Packaging

RG5401



Features Include

- Strong Lightweight Aluminium Construction
- High User Capacity 159kg
- Ergonomically Designed Seat for Maximum Comfort
- Adjustable Aluminium Legs with Sure Lock Adjustment & Flexible PVC Leg Tips
- Easy Assembly (**No Tools Required**)
- Supplied in KD Form for Easy Transport
- Ideal for Internet Sales
- Easy Clean Design
- Packaged for Display

Product Dimensions

Frame Width	485mm
Seat Width	435mm
Seat Depth	385mm
Seat Handle Height	160mm - 180mm
Seat Height	465mm to 575mm
Maximum User Weight	159kg
Carton Dimensions	620 x 500 x 20mm 4.0kg

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@shoprider.com.au