

Folding Over Seat Toilet Riser



Features

- Lightweight Aluminium Construction
- Pearl White Paint Finish
- Easy Folding for Storage
- Easy Clean Hygienic Seat and Lid
- Includes Splash Guard & Bowl
- Multi Slot Legs for Uneven Floors
- High Grade Rubber Feet
- Rust Free Fittings

Dimensions

- Seat Width - 480mm
- Seat Height - 460mm - 620mm
- Max User Capacity - 100kg
- Weight - 3.5kg

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@shoprider.com.au