

Shower Boards

RG9201 / RG9202



Features

- Lightweight polyurethane material
- Textured non slip surface
- Easy clean surface
- Drain holes in seat
- Integrated soap holder
- Non slip adjustable feet
- Anti slip rubber pads on the base

Specifications

RG9201

- Max User Capacity - 150kg
- Width Adjustment - 450 to 710mm
- Seat Width - 387mm

RG9202

- Max User Capacity - 150kg
- Width Adjustment - 400 to 630mm
- Seat Width - 350mm

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@shoprider.com.au