

# Excel

# !CLICK-&-GO Lite™



Innovation in wheelchair power support systems

Compact powerpack

Simply to fix under the wheelchair

Easy to operate



The Excel Click & Go Lite is equipped with two drive wheels and programmable microprocessor-driven electronics. The unit can be fitted to most wheelchairs with seat widths between 40 and 55 cm. Easy operating at the push handles. The Click & Go Lite is light enough to take anywhere you go, and when you do not need it, the pack is easily lifted off the floor.



Small enough to take everywhere you go.



Always smart: easy to adjust speed forwards and backwards drive.



Easy to fit

‘Smart mobility solutions’

Maximum speed:	4 mph
Range:	8 km
Motor:	12V, 90W-280W
Battery:	12V 20Ah
Charger:	3A/12V

Fits wheelchairs seat width with:	40 to 55 cm
-----------------------------------	-------------

Overall width:	25 cm
Overall length:	44 cm
Overall height:	29½ cm
Drive Wheels:	200x50 mm

Weight (complete):	14 kg
Weight (without batteries):	7,7 kg
Weight (battery):	6,3 kg

The Click & Go comes standard with:  
- All fittings needed to get you going  
- Battery and charger  
- Hand control unit



We hold a copyright protection (patent) and consequently reserve the right to amend design and specifications without prior notice all errors and omissions excepted. All rights reserved. Copyright 2006 Van Os Medical B.V.

[www.vanosmedical.com](http://www.vanosmedical.com)



Van Os Medical UK Ltd

Excel House  
Escrick Business Park  
Escrick, York  
YO19 6FD

Tel: 01904 720170  
Fax: 01904 720398

[sales@vanosmedical.com](mailto:sales@vanosmedical.com)