

Shoprider

889 GOLF

ASK FOR THE BRAND SHOPRIDER



889 GOLF

SHOPRIDER[®]
mobility
www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER[®] *mobility* **889 GOLF**

www.shoprider.com.au

The SHOPRIDER™ **889 GOLF** has been designed to take on the challenges of the Golf course while providing the reliability and comfort of the SHOPRIDER range. The **889 GOLF** comes with wide tyres 330 x 120, 2.0hp High torque motor, High back seat and all the features you expect from the SHOPRIDER range, Fully rotating lockable seat, suspension, adjustable steering tiller and multi position fold down arm rests.

The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS	889 GOLF
Length (*Including Golf Bracket)	1430mm (*1460mm)
Width	750mm
Height	1300mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
*Maximum Load Single user	225Kgs
Charger	8Amp
Motor Size	2.0Hp
Maximum Speed	10Kph
Battery capacity (Minimum Amp)	2 x 55 Amp
Maximum range#	up to 45km
Wheel Size	330mm X 120mm
Turning Radius	1800 mm
Weight (with battery)	146kg

FEATURES:

- High Back Seat
- T-Bar Steering Tiller
- Adjustable Steering Tiller
- Adjustable width arm rests
- Sliding Rotating Seat
- Variable Seat Height
- Battery Power indicator
- Double Head Lights
- Rear Lights & Indicators
- Hazard Lights
- Front & Rear Bump Pads
- High & Low Speed
- * Golf Bracket can be fitted



Comprehensive Warranty


Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

WARNING: Do not exceed 10kph on public roads or paths

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



SHOPRIDER[®] *mobility*

www.shoprider.com.au

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au