

Shoprider AllRounder

ASK FOR THE BRAND SHOPRIDER



SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER® AllRounder

www.shoprider.com.au *mobility*

The SHOPRIDER™ **AllRounder** is the new large scooter from SHOPRIDER and it is loaded with great features. The **AllRounder** has all the standard SHOPRIDER features 330mm tyres for safety & comfort, fully rotating lockable seat, adjustable steering tiller and fold down arm rests. The reliability of these impressive scooters by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS

	AllRounder
Length	1300mm
Width	640mm
Height	1200mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
Maximum Load Single User*	159Kgs
Charger	8Amp
Motor Size	1.5Hp
Maximum Speed	12Kph
Battery Capacity (Minimum Amp)	40Amp
Maximum Range #	up to 40km
Wheel Size	330mm
Turning Radius	1440mm
Weight (with battery)	99kg

FEATURES

- Large Captain's Seat
- Adjustable Steering Tiller
- Adjustable Width Arm Rests
- Sliding rotating Seat
- Fold Arm Rests
- Variable Seat height
- Front Shopping Basket
- Horn
- Led Battery Level indicator
- Front Chrome Bumper
- Front & Rear Suspension
- Front & Rear Lights
- Hazard Lights /indicators
- High / Low Speed switch
- Anti Tips Wheels
- Comprehensive Warranty**

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)


WARNING: Do not exceed 10kph on public roads or paths.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

AllRounder



SHOPRIDER®
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au