

Foam Cushion with Gel Insert



This combination foam and gel cushion is suited to people at low to medium risk of pressure injury, and helps to relieve pain from prolonged sitting.

Made from a gel bladder encased inside a resilient polyurethane foam cushion, this design helps the cushion to evenly disperse the weight of a seated person over a larger area, which minimises pressure points on the buttocks or thighs.

The water based gel not only helps reduce direct pressure on the skin that might lead to pressure injury, but also provides a cooling effect to reduce heat and perspiration build up. The high density foam provides stability plus great comfort for prolonged sitting.

The outer cover is low shear to prevent damage to skin as a result of friction – the damage that can be caused to delicate skin when it is dragged across another surface.

It also has a non-slip base to help stop the cushion from sliding off the chair. Two positioning handles at either side of the cushion make it easy for the user to reposition the cushion.

CODE	BL9375	BL9382	BL9399
HEIGHT	50 mm	50 mm	50 mm
WIDTH	400 mm	450 mm	450 mm
DEPTH	450 mm	450 mm	500 mm
CAPACITY	110 kg		
WARRANTY	1 year		

