

# Silicone Fibre Cushion

Cotton/polyester breathable cover provides maximum cool comfort



Silicone coating lubricates each fibre to reduce friction and shear force while preventing moisture penetration or retention



This seat cushion is constructed from silicone coated cylindrical fibres that retain their loft and promote airflow to the skin's surface for superior patient support and comfort.

The silicone coating lubricates each fibre to reduce friction and shear force while preventing moisture penetration or retention, making it an ideal choice for geriatric or palliative care.

Covered with a breathable cotton/polyester fabric, the cushion is fully machine washable with standard detergent.

- Excellent pressure relief and comfort on any chair or seat surface
- Fluid resistant silicone fibres provides extra support and comfort
- Ideal for geriatric and palliative care
- Flame Retardant 16 CFR Part 1632
- Fully machine washable up to 70°

CODE	BL7302	BL7319	BL7326	BL9979
SIZE	16 x 18"	18 x 16"	20 x 16"	20 x 17"
HEIGHT	80 mm	80 mm	80 mm	80 mm
DEPTH	460 mm	410 mm	410 mm	430 mm
WIDTH	410 mm	460 mm	510 mm	510 mm
CAPACITY	110 kg			
CONFORMS TO	ARTG176835			
WARRANTY	1 year			

