

The BetterLiving alternating pressure surfaces is designed for the prevention and treatment of skin breakdown and pressure injuries.

		Mattress Replacement	Mattress Overlay	
General	CODE	APMBL-R01	APMBL-L01	
	CAPACITY (SAFE WORKING LOAD)	40-180 kg	40-150 kg	
	SYSTEM ALTERNATION	1 in 2		
	WARRANTY	1 year		
	COMPLIANCE	IEC 60601-1, IEC60601-1-2, IEC60601-1-11		
	ARTG	289458		
	PRESSURE INJURY RISK	High to very high risk	At risk to high risk	
	CARTON DIMENSIONS (L x W x H) in mm	470 x 470 x 430 mm		
	DIMENSIONS (L x W x H) in mm	2000 x 880 x 200 mm	2000 x 880 x 130 mm	
	MATTRESS WEIGHT	8 kg	4.5 kg	
Support Surface	CPR TYPE	QuickTwist CPR		
	CELLS	NUMBER	18	
		STATIC HEAD CELLS / ZONE	3	
		CELL-IN-CELL	Yes	No
		CELL MATERIAL	TPU-laminated nylon	
		MATTRESS AIR SUPPLY LINE	Heavy duty PVC piping	
		MATERIAL / TYPE	Polyurethane-laminated nylon	
		FIRE RETARDANCY	BS7175 Crib 5	
	TOP COVER	MAXIMUM WASH TEMPERATURE	95° C	
		ATTACHMENT	Dual zippers	
		SEAM CONSTRUCTION	All external seams are welded for infection control	
	BASE	MATERIAL / TYPE	PVC-laminated Polyester	
		FIXING STRAPS	Corner fastening straps	Head/foot elastic straps

Optional  
Control  
Unit

<b>DIMENSIONS (L x W x H) in mm</b>		120 x 265 x 195 mm
<b>WEIGHT</b>		2 kg
<b>POWER</b>		100-240V 50/60Hz
<b>COMPRESSORS</b>		Single
<b>OPERATION</b>		Full digital control
<b>CYCLE TIME</b>		12 minutes
<b>NOISE LEVEL, WITHIN 2 METRES</b>		32 dB(A) ±1
<b>MIN / MAX PRESSURE</b>		20-60mm HG +/- 6mm HG
<b>ELECTRICAL PROTECTION TYPE</b>		Class II Type BF
<b>OPERATING ENVIRONMENT</b>		Operation 15°C to 35°C Storage 5°C to 60°C
<b>PRESSURE SETTING</b>	<b>AUTOMATIC</b>	Yes
	<b>MANUAL</b>	5 levels
<b>ALARMS</b>		Pressure Loss Power Failure



THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.