

Alternating Mattress / Overlay



BetterLiving alternating pressure surfaces are designed to offer a cost-effective solution for the prevention and management of pressure related injuries.

Harnessing the latest technical features for improved clinical care and ease of use, BetterLiving offers both a mattress replacement and a mattress overlay option. Ideally suited to home care or aged care environments, these products can be used for people at risk to high risk of pressure injury.

ITEM	Mattress Replacement	Mattress Overlay
CODE	APMBL-R01	APMBL-L01
FEATURES	<ul style="list-style-type: none"> Essential features for pressure injury prevention and treatment Automatic pressure setting Replace existing mattress for superior therapy Cell in cell for enhanced user safety and comfort Easy to operate and maintain 	<ul style="list-style-type: none"> Functional overlay for pressure prevention Automatic pressure setting Used over existing mattress for basic support* Easy to operate and maintain
CAPACITY	40 - 180 kg	40 - 150 kg
PRESSURE INJURY RISK	High to very high risk	At risk to high risk
COMPLIANCE	IEC 60601-1, IEC60601-1-2, IEC60601-1-11	
ARTG	289458	
WARRANTY	1 year	

* Recommended for use with a Safety Underlay; must be used with a mattress

BetterLiving is a brand of cost effective yet clinically effective products designed for the Australian healthcare market. Our extensive product range includes beds and accessories, bathroom and hygiene, seating and positioning, mobility equipment and a range of aids for daily living.

As specialists in the supply of pressure care and manual transfer equipment solutions, Novis is pleased to launch BetterLiving's first range of affordable alternating pressure surfaces.

Our focus is to offer the most efficient products at the most cost effective price – without compromise to product integrity or clinical patient care.

Alternating Mattress/ Overlay



Why choose BetterLiving?

Safety

No more concerns over correct pressure level settings with our AutoCair feature – automatic weight sensing technology that delivers fail safe patient comfort levels irrespective of their position.

With the AutoCair feature active, the system adjusts to the appropriate pressure setting based on the patient's weight, reducing the need for carer intervention. Manual adjustments can still be made for individualised therapy or comfort.

Simplicity

Offering all the features you need, the BetterLiving pressure surfaces are designed for basic ease of use and minimal maintenance. Easy to set up, simple to operate and straightforward to clean and maintain, these products require minimal technical training or support.

Quality

Our BetterLiving pressure surfaces conform to all relevant Australian standards. Made using the best quality materials, every product must pass our stringent Quality Control inspection and regular testing protocols.

