

Wideback Backrest Support



Fabric: open weave,
non-allergenic
brushed nylon


Core: three grades
of high-resilient
urethane foam

Frame: impact resistant
polycarbonate
biomechanically moulded



Specially designed for
broader backs or larger frames

Removable, adjustable
lumbar pad



An award-winning portable backrest that transforms ordinary chairs into ergonomically correct seating. This wider support is specifically designed for users with a broad back or larger build.

Unlike ordinary "L" shaped chairs, the "S" shape of this backrest moulds the spine into a more anatomically correct position, helping to enhance overall posture by supporting proper spinal alignment.

- designed for users with a broad back and/or generally larger frame (suits male size 44+ and female size 18+)
- 50 mm wider and 30 mm taller than the Lowback Backrest Support
- lightweight and portable; includes elastic positioning straps and handle
- unique polycarbonate frame with polyurethane foam gently supports natural spine contours
- removable and adjustable lumbar pad for additional, customisable back support
- hypoallergenic brushed nylon cover can be removed and sponge washed
- helps prevent fatigue, provide relief from pain aggravated by poor posture and ease pain caused by spinal disc problems and sciatica

CODE OBUSWB

COLOUR Black

WIDTH 500 mm

HEIGHT 560 mm

DEPTH 75 mm

MATERIALS Nylon, Foam,
Polycarbonate

WEIGHT 2 kg

CAPACITY 200 kg

WARRANTY 1 Year

ENDORSED BY

