

Introducing ProCair  
by Novis – alternating  
solutions for the prevention and  
treatment of pressure injuries and other  
preventable complications of immobility.

		Standard size	
General	<b>CODE</b>	APMPC-R03	
	<b>CAPACITY (SAFE WORKING LOAD)</b>	30-200 kg	
	<b>SYSTEM ALTERNATION</b>	1 in 2	
	<b>WARRANTY</b>	3 year	
	<b>CARRY BAG</b>	Durable black carry bag	
	<b>PRESSURE INJURY RISK</b>	High to very high risk	
	<b>COMPLIANCE</b>	IEC 60601-1, IEC60601-1-2, IEC60601-1-11	
	<b>ARTG</b>	289458	
	<b>CARTON DIMENSIONS (L x W x H) in mm</b>	1000 x 670 x 270 mm	
	<b>DIMENSIONS (L x W x H) in mm</b>	2000 x 900 x 180 mm	
Support Surface	<b>MATTRESS WEIGHT</b>	12 kg	
	<b>CPR TYPE</b>	QuickTwist CPR at head end	
	<b>CELLS</b>	<b>NUMBER</b>	21
		<b>STATIC HEAD CELLS / ZONE</b>	3
		<b>NARROW FOOT CELLS / ZONE</b>	4
		<b>BOTTOMING-OUT PROTECTION</b>	In-built foam insert and foam sheet
	<b>CELL MATERIAL</b>	TPU coated nylon	
	<b>TOP COVER</b>	<b>MATERIAL / TYPE</b>	Dark blue, polyurethane-laminated nylon – water resistant, vapour permeable multi-stretch
		<b>FIRE RETARDANCY RATE</b>	BS7175 Crib 5
		<b>MAXIMUM WASH TEMPERATURE</b>	95° C
		<b>ATTACHEMENT</b>	Dual zippers
		<b>SEAM CONSTRUCTION</b>	RF welded on all external seams
		<b>WATERFALL SKIRT</b>	Double layered
		<b>MATERIAL / TYPE</b>	Black PVC laminated polyester
	<b>BASE</b>	<b>HINGE</b>	No
<b>FIXING STRAPS</b>		Webbing at 4 corners, 4 at middle	
Therapeutic Modes	<b>DYNAMIC MODE</b>	Alternating air cells allowing complete pressure offload and optimal therapeutic benefits. Foam sheet under air cells provides bottoming out protection.	
	<b>HYBRID MODE</b>	Foam sheet placed over air cells to allow immersion and greater comfort	

# NOVIS

## ProCair Prime

### TECHNICAL SPECIFICATIONS

Built-in  
Control Unit

<b>POWER</b>	100 - 240V 50/60Hz	
<b>COMPRESSORS</b>	Single	
<b>OPERATION</b>	Full digital control	
<b>CYCLE TIME</b>	12 minutes	
<b>NOISE LEVEL, WITHIN 2 METRES</b>	32 db(A) ±1	
<b>MIN / MAX PRESSURE</b>	20-60mm HG +/- 6mm HG	
<b>ELECTRICAL PROTECTION TYPE</b>	Class II Type BF	
<b>OPERATING ENVIRONMENT</b>	Operation 15°C to 35°C Storage 5°C to 60°C	
<b>PRESSURE SETTING</b>	<b>AUTOMATIC</b>	Yes - cyclical, self-monitors and adjusts every cycle
	<b>MANUAL</b>	5 levels
<b>ALARMS</b>		Low Pressure
		Power Failure



THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL



Pressure care and patient handling specialists

novis.com.au / 1300 738 885