

Rollator - W463

Features

Lightweight aluminium frame

Durable 8" puncture free foam tyres

Easy to fold for storage & transportation

Padded seat with removable folding backrest



colour available - Blue or Red

Technical Specs

| Seat Width | 35cm |
|------------------------|---------|
| Seat Depth | 33cm |
| Seat Height | 53 cm |
| Width at base | 69cm |
| Depth at base | 64cm |
| Height | 83-96cm |
| Weight Capacity | 114kg |
| Total weight | 7kg |
| | |